



PRODUCT DATA SHEET

Page 1 / 2

PBS Endotoxin-free (sterile) [For Nano-formulated Drug Analysis]

Cat. No.: IAX-900-001DC **Lot. No.:**

Synonyms	Phosphate-buffered saline (PBS) Endotoxin-free (sterile)		
Formulation	10ml (L010), PBS in double-distilled water (ddWater), without calcium and magnesium, particle-test		
Endotoxin-free	Bacterial Endotoxin Test (kinetic turbidimetric LAL method) according to Ph. Eur. 9. Passed according to specification: Endotoxin-free: <0.00002 [EU/mg]/<0.02 [EU/ml].		
Sterility	Filter method: according to Ph. Eur. 9. Passed according to specification: No growth in Thioglycolate medium at 30-35°C after 14 days. No growth in Soybean Casein Digest Broth (TSB) at 20-25°C after 14 days.		
pH Value	7.4 (± 0.2)		
Composition	Sodium chloride: 137mM, Potassium chloride: 2.7mM, Potassium phosphate monobasic: 2mM, and Sodium phosphate dibasic: 10mM in ddWater. Does not contain calcium or magnesium.		
Application	PBS Endotoxin-free (sterile) is a balanced salt solution used for a wide range of applications (cell culture, washing cells, transporting cells or tissue, diluting cells for counting, and preparing other laboratory reagents). This PBS product is without calcium and magnesium for the purpose of rinsing cells before treating them with proteases for dissociation.		
Handling	Keep sterile. Ready-to-use, sterile solution. Avoid plastic containers, or any source of MNP.		
Shipping	Ambient		
Storage	2-8°C		
Stability	 Upon receipt, store product at 2-8°C. Do not freeze. In its unopened original vial, the product is stable for at least 24 months when stored at 2-8°C. Once the vial is opened, or if aliquoted into sterile vials under sterile conditions, the product remains stable for an additional 12 months at 2-8°C. 		
MSDS	Available on request		

Document No.: | AX-900-001DC | **Version:** | 1.2 | **Issue Date:** | 13/11/2024

DISCLAIMER: THIS PRODUCT IS NOT INTENDED OR APPROVED FOR HUMAN, DIAGNOSTICS OR VETERINARY USE. USE OF THIS PRODUCT FOR HUMAN OR ANIMAL TESTING MAY BE EXTREMELY HAZARDOUS AND MAY RESULT IN DISEASE, SEVERE INJURY, OR DEATH. THIS PRODUCT IS FOR RESEARCH USE ONLY (RUO).





PRODUCT DATA SHEET

Page 2/2

PBS Endotoxin-free (sterile) [For Nano-formulated Drug Analysis]

Cat. No.: IAX-900-001DC **Lot. No.:**

Endotoxin-free and Sterile I	Buffers and Rela	ted Products
------------------------------	------------------	--------------

IAX-900-001	PBS Endotoxin-free (sterile)
IAX-900-001DC	PBS Endotoxin-free (sterile) [For Nano-formulated Drug Analysis]
IAX-900-002	ddWater Endotoxin-free (sterile)
IAX-900-002DC	ddWater Endotoxin-free (sterile) [For Nano-formulated Drug Analysis]
IAX-900-003	Physiological Saline [Sodium Chloride 0.9% Endotoxin-free] (sterile)
IAX-900-003DC	Physiological Saline [Sodium Chloride 0.9% Endotoxin-free] (sterile) [For Nano-formulated Drug Analysis]
IAX-900-004	PBS with EDTA Endotoxin-free (sterile)
IAX-900-005	TRIS with EDTA [TE Buffer] (100x) Endotoxin-free (sterile)
IAX-900-006	EDTA (400mM) Endotoxin-free (sterile)
IAX-900-007	HEPES Buffer (500mM) Endotoxin-free (sterile)
IAX-900-008	DNA Loading Buffer with TRIS and EDTA (6x) (Blue)
IAX-900-009	HEPES Buffer (50mM) with NaCl [Sodium Chloride] (150mM) Endotoxin-free (sterile)
IAX-900-010	NaCl [Sodium Chloride] (1.5M) Endotoxin-free (sterile)
IAX-900-011	TRIS Buffer (I.5M) Endotoxin-free (sterile)
IAX-900-012	TRIS Buffer (30mM) with NaCl [Sodium Chloride] (150mM) Endotoxin-free (sterile)
IAX-900-013	PBS with Magnesium and Calcium Endotoxin-free (sterile)
IAX-900-014	ddWater with 0.9% Benzyl Alcohol [Bacteriostatic Water] Endotoxin-free (sterile)

Document No.: | AX-900-001DC | **Version:** | 1.2 | **Issue Date:** | 13/11/2024